

Application Note

**Computing Forming-Limit Curves
(FLC ISO 12004) in VIC-3D**

VIC-3D 11

2026

VIC-3D 11

Computing Forming-Limit Curves (FLC ISO 12004) in VIC-3D 11

This application note describes the recommended workflow for computing forming-limit strain pairs (ϵ_1 , ϵ_2) from DIC measurements using two Correlated Solutions extensions for VIC-3D:

FLC ISO 12004

This per-sample analysis implements the position-dependent method of **DIN EN ISO 12004-2:2021** for a single fractured specimen and returns one (ϵ_1 , ϵ_2) limit-strain pair per section line.

FLC ISO 12004 Annex F

This aggregator combines the per-sample CSV outputs into the consolidated table prescribed by **Annex F** of the standard, with means per-sample and per-geometry.

Please note: This note assumes you are already following the standard for specimen preparation, speckle application, and test execution for your laboratory. It does not provide a step-by-step walkthrough of DIC (test set-up, calibration, image acquisition, normal post-processing, etc.). For a range of resources focused on DIC acquisition and analysis, please visit the Correlated Solutions Support Page at www.correlatedsolutions.com/support or by clicking the QR code to the right.



Before You Begin

Each sample project must contain:

- A completed correlation through fracture
- The working directory named in the pattern 'X.Y.Z' (e.g. '1.20.1'), where the middle field is the **sample geometry** and the last is the **sample number**. Use numeric digits in all three fields — non-numeric values will break the Annex F aggregation.
- A single inspector line drawn on the analysis frame (see *Step 2*).

The screenshot shows a file explorer view of a project directory. The root directory is 'project_root/'. Inside, there are four sub-directories representing different sample projects, each named with a pattern 'X.Y.Z'. The sub-directories are: '1.20.1/' (geometry '20', sample '1'), '1.20.2/' (geometry '20', sample '2'), '1.40.1/' (geometry '40', sample '1'), and '1.40.2/' (geometry '40', sample '2'). Each sub-directory contains a 'VIC' folder (VIC-3D project files) and a 'flc_results.csv' file. At the root level, there is a 'flc_annex_F.csv' file, which is the aggregated output.

Directory naming pattern

1 . 20 . 1 /

sample number
sample geometry

- All three fields must be numeric digits.
- Other names break Annex F grouping.


Legend

- Parent directory — run the Annex F extension from any sibling.
- Sample project — one per replicate, named X.Y.Z.
- 📄 Per-sample result — written by FLC ISO 12004.
- 📄 Aggregated output — written by FLC ISO 12004 Annex F.

Step 1 - Load the Analysis Frame

The position-dependent method assumes the localized neck is fully developed and the crack has **not yet** propagated through the gauge region.


Step backward through the test until you find the **last frame before crack-through**. Note the frame number, this will be the data file selected within the extension.



Step 2 - Draw the Inspector Line


Draw one inspector line across the cracked region. It should be:

- Perpendicular to the crack (along the direction of maximum ϵ_1).
- At least 40 mm long. The extension rejects lines shorter than 20 mm and §7.2.2 of the standard requires a 40 mm section length. Draw with at least 20 mm of usable data on each side of the visible crack.





Step 3 - Run the FLC ISO 12004 Extension

- The FLC ISO 12004 Extension default settings are set according to the standard.
- Options are well documented within the extensions Help page.
- The only inputs required will be the 'File' that was noted previously as the last frame before crack-through, and the inspector item drawn across the crack in Step 2.
- Complete this step for all of the samples in the series.



Once this is complete, **Section lines (right)** should be visible in the data viewer, centered on the crack and perpendicular to it.



Additionally, Line Extractions are added to the project. These contain the ϵ_1 and ϵ_2 profiles for each section line. Each ϵ_1 profile should show a single sharp peak as shown in the example (left).

In addition to these new values within the project there will also be two CSV files written to the working directory:

- **flc_profile_data.csv** — full ϵ_1 and ϵ_2 profiles for every section line
- **flc_results.csv** — one (ϵ_1 , ϵ_2) limit-strain pair per successfully analyzed section line.

Verify that the **Sample geometry** and **Sample number** columns contain the expected values from your directory name. If they read 'NA', the directory does not match the 'X.Y.Z' pattern. If this is the case, rename it and rerun.

Step 4 - Aggregating Multiple Samples - FLC ISO 12004 Annex F

Upon execution of the FLC ISO 12004 extension on every sample, use the Annex F extension to combine the results.

All sample projects must sit as siblings under a common parent directory. Each directory name must follow the 'X.Y.Z' pattern with **numeric digits** in all three fields — the middle field groups by geometry; the last groups duplicates within a geometry.

Directory naming pattern

1.X.Y.Z

- All three fields must be numeric digits.
- Other names break Annex F grouping.

Legend

- Parent directory — run the Annex F extension from one sibling.
- Sample project — one per replicate, named X.Y.Z.
- 📄 Per-sample result — written by FLC ISO 12004.
- 📄 Aggregated output — written by FLC ISO 12004 Annex F.

Open any of the sample projects (or a dedicated summary project in the same parent), then **Run** the extension *Extensions* → *Forming limit curve* → *FLC ISO 12004 Annex F*.

Executing the extension triggers the following actions:

- Scan of all sibling directories in the current project for **flc_results.csv**.
- Concatenation of all rows into a single table.
- Computation of the **mean ϵ_1 and ϵ_2 per-sample** (grouped by sample geometry + sample number).
- Computation of the **mean ϵ_1 and ϵ_2 per-geometry** (grouped by sample geometry).
- Write the aggregated CSV to the parent directory.

The per-sample and per-geometry mean columns appear on the last row of each sample / each geometry, matching the Annex F table format.

ϵ_1	ϵ_2	ϵ_1 - Mean sample	ϵ_2 - Mean sample	ϵ_1 - Mean geometry	ϵ_2 - Mean geometry	Sample geometry	Sample number	Section line
0.494407922	-0.225146781						1	1
0.506876114	-0.223545249						1	1
0.509517531	-0.218689753	0.502319698	-0.217805425				1	1
0.513156183	-0.222643733						1	2
0.510168945	-0.221764868						1	2
0.505709023	-0.218790681	0.508633404	-0.224630386				1	2
0.504596419	-0.227966238						1	3
0.514275835	-0.218908373						1	3
0.508939729	-0.219885998	0.500240587	-0.221121355	0.507223338	-0.223956257		1	3
0.50700477	-0.196351904						2	1
0.50974718	-0.190808476						2	1
0.515660284	-0.196738648	0.517861947	-0.19449056				2	1
0.507910866	-0.192325528						2	2
0.500300642	-0.187504554						2	2
0.509630786	-0.188902173	0.511095606	-0.188715061				2	2
0.505007439	-0.194574429						2	3
0.500668963	-0.194513781						2	3
0.505139678	-0.196753829	0.505402118	-0.193219563	0.498657334	-0.189648007		2	3
0.49219864	-0.149559934						3	1
0.477852984	-0.151189535						3	1
0.495020844	-0.152946243	0.483870237	-0.149331878				3	1
0.489586922	-0.154504271						3	2
0.492980555	-0.153637104						3	2
0.4848902	-0.155433138	0.483710037	-0.154981313				3	2
0.479746162	-0.150662225						3	3
0.478401955	-0.149486863						3	3
0.495532284	-0.154415315	0.483566701	-0.149531801	0.476313346	-0.148869363		3	3
0.468239144	-0.109827418						4	1
0.471797648	-0.110587062						4	1
0.47565772	-0.11075452	0.480388754	-0.113552745				4	1
0.467065382	-0.110620216						4	2
0.47853246	-0.10831976						4	2
0.466576395	-0.109573597	0.466987543	-0.108948812				4	2
0.470888475	-0.107964886						4	3
0.470888413	-0.109463485						4	3
0.46474947	-0.107399829	0.464888696	-0.108539252	0.473714219	-0.107738746		4	3
0.473107947	0.44404994						8	1
0.470867027	0.45351939	0.466224424	0.446863907				8	1
0.479000504	0.453032846						8	3
0.477430641	0.462096254						8	3
0.465308673	0.460181149	0.481602157	0.451110939	0.475176919	0.44949313		8	3

Support

If you have any questions about this document or any other questions, comments, or concerns about our software, please contact us at support@correlatedsolutions.com, or visit our website at correlatedsolutions.com/support.

The powerful extensions feature is available exclusively in VIC-3D 11.2. Contact one of our Sales Engineers at sales@correlatedsolutions.com for information on updating your software.